


LABTECH INTERNATIONAL

| | |
|---------------------------|---|
| Item Name : | Centrifuge Tubes, Conical Bottom Without Graduation |
| Item Code : | 7680 |
| Item Description : | <ul style="list-style-type: none"> • Made from ASTM E-438 Type-1 Class A, Boro 3.3 glass • Uniform wall Thickness • Strudy Lip |
| Item Image : |  |

| Part No. | Capacity (ml) | Height x O.D. (mm) | PACK Qty. |
|----------|---------------|--------------------|-----------|
| 7680-5 | 5 | 100 x 13 | 100 |
| 7680-10 | 10 | 110 x 15 | 100 |
| 7680-15 | 15 | 120 x 17 | 100 |
| 7680-25 | 25 | 125 x 22 | 100 |
| 7680-50 | 50 | 125 x 29 | 50 |

11-A, Herman City Road, Khudda Khurd, Ambala Cantt, Haryana 133006

+91-9802013010, +91-9896317240

info@labtechexport.com, <http://www.labtechexport.com>