



## LABTECH INTERNATIONAL

<b>Item Name :</b>	Centrifuge Tubes Box PP
<b>Item Code :</b>	7730
<b>Item Description :</b>	

- Made of Polypropylene
- These boxes provide substitute for conventional cardboard or foam rack
- These boxes are autoclavable & can withstand Temperature ranging from -90°C to +120°C

<b>Item Image :</b>	
---------------------	---

Part No.	Capacity (ml)	Pack Qty	
7730-30	30	4	
7730-50	50	4	

11-A, Herman City Road, Khudda Khurd, Ambala Cantt, Haryana 133006

+91-9802013010, +91-9896317240

info@labtechexport.com, <http://www.labtechexport.com>