



LABTECH INTERNATIONAL

Item Name :	ADAPTERS, STOPCOCK, SOCKET
Item Code :	2560
Item Description :	<ul style="list-style-type: none">•€With PTFE needle valve stopcock providing fine flow control all positive closure•€Detachable plastic screw thread connectors allows hose etc.•€To be fitted easily and safely
Item Image :	

PART No.	NS FEMALE SIZE	STOPCOCK	Pack(Qty)	
2560-14F	14/23	3.0	10	
2560-19F	19/26	3.0	10	
2560-24F	24/29	3.0	10	

11-A, Herman City Road, Khudda Khurd, Ambala Cantt, Haryana 133006

+91-9802013010, +91-9896317240

info@labtechexport.com, <http://www.labtechexport.com>