

## LABTECH INTERNATIONAL

Item Name :	SCREWTHREAD, GL
Item Code :	1410
Item Description :	

- Externally threaded glass connectors.
- For use as replacement on apparatus using the GL thread or when designing items where the external thread is preferred.

Item Image :	
--------------	---

Part No.	Thread size (G.L No.)	Tubing O.D.Xwall,mm	Pack Qty.
1410-14GL	GL-14	12x1.5	100
1410-18GL	GL-18	16x1.8	100
1410-25GL	GL-25	22x1.8	100
1410-32GL	GL-32	28x2.0	50
1410-45GL	GL-45	40x2.3	10
1410-120GL	GL-120	120x5.0	10

11-A, Herman City Road, Khudda Khurd, Ambala Cantt, Haryana 133006

+91-9802013010, +91-9896317240

info@labtechexport.com, <http://www.labtechexport.com>