



## LABTECH INTERNATIONAL

<b>Item Name :</b>	CLEAR GLASS VOLUMETRIC FLASK WITH SOLID GLASS STOPPER, CLASS A, SERIALIZED WITH INDIVIDUAL CERTIFICATE ASTM E-288
<b>Item Code :</b>	2740
<b>Item Description :</b>	

**Volumetric flasks are made to the highest standards to qualify for Class A certification.  
Flasks are individually serialized.  
Made from Borosilicate 3.3 expansion glass.  
Flasks are serialized and supplied with glass PENNY HEAD STOPPER.**



<b>Item Image :</b>	
---------------------	--

Part No.	Capacity (ml)	±Tol.(ml)	Stopper	Pack QTY	
2740-5	5	0.02	9	10	
2740-10	10	0.02	9	10	
2740-20	20	0.03	9	10	
2740-25	25	0.03	9	10	
2740-40	40	0.04	9	10	
2740-50	50	0.05	9	10	

2740-100	100	0.08	13	10	
2740-200	200	0.10	13	10	
2740-250	250	0.12	16	10	
2740-500	500	0.20	19	10	
2740-1000	1000	0.30	22	10	
2740-2000	2000	0.50	27	10	

11-A, Herman City Road, Khudda Khurd, Ambala Cantt, Haryana 133006

+91-9802013010, +91-9896317240

info@labtechexport.com, <http://www.labtechexport.com>